Experimental Technique for Radiative-Process-Resolved X-Ray Absorption Spectroscopy at the Inner-Shell Excitation Thresholds

Y. Muramatsu¹, and R. C. C. Perera²

¹Japan Atomic Energy Research Institute, Sayo, Hyogo 679-5148, Japan

²Center for X-Ray Optics, Lawrence Berkeley National Laboratory, Berkeley, California 94720, USA

Partial-fluorescence-yield (PFY) x-ray absorption measurements, using the optimized window widths of position sensitive detectors in wave-length dispersive x-ray spectrometers, have been applied for radiative-process-resolved (RPR) x-ray absorption spectroscopy [1]. For such the PFY absorption measurements, we can selectively extract x-ray emission signals from specified radiation processes by optimizing the window width on the PSD of a grating x-ray spectrometer installed in BL-8.0.1. Figure 1 shows the top view of the optical arrangement for the measurements. In this arrangement, the high-energy and low-energy window widths on the PSD can be varied by adjustments of the high- and low-energy masks. To measure the PFY absorption spectra, we scanned the incident beam energy monitoring the x-ray intensities detected on the PSD whose window widths were optimized by its high- and low-energy masks. We measured PFY-absorption spectra of graphite and diamond at the C K threshold and of h-BN and c-BN at the B K threshold. From these measurements, resonant elastic x-ray scattering was observed in graphite and h-BN on their PFY absorption spectra, and excitonic x-ray scattering We also measured the PFY-absorption spectrum of was observed in diamond and c-BN [1]. NiO at the Ni L threshold, resolving the L α and L β x-ray emission processes, as shown in Figure 2. These results show that PFY absorption measurements for RPR x-ray absorption spectroscopy can provide the information about the electronic structures and the radiative-decay process in inner-shell excitation.




Figure 1. Top view of the optical arrangement for the PFY absorption measurements in BL-8.0.1.




Figure 2. Ni L x-ray emission spectrum (upper panel) and PFY and total-electron-yield (TEY) x-ray absorption spectra (lower panel) of NiO. The PFY x-ray absorption spectra were taken by monitoring the x-ray intensity measured by the PSD with (A) wide, (B) low-energy, and (C) high-energy windows. Spectrum obtained after subtracting (B) from (A), denoted by (A) - (B) is also shown.

[1] Y. Muramatsu et al., J. Synchrotron Radiation (in press).

This work was supported by the Japan Atomic Energy Research Institute.

Principal investigator: Yasuji Muramatsu, Japan Atomic Energy Research Institute. Email: murama@spring8.or.jp. Telephone: +81-791-58-2601.